

【英文名称】: TUMOUR MARKER CONTROL - LEVEL 3 (TMR CONTROL 3)

【中文名称】: 肿瘤标记物质控品（水平 3）

【货号】 (CAT. NUMBER): TU5003

【批号】 (LOT NUMBER): 409TU

【包装规格】 (SIZE): 3 x 2 ml

【到期】 (EXPIRY): 2027-04-28

本产品的准备和使用详见试剂盒内说明书。

Full instructions for the use and preparation of this product are included with the kit.

质控品靶值和范围如下所示。

Target values and ranges are listed below.

ZH: 23 Apr '25 me

Method 所有仪器均值		TUMOUR MARKER CONTROL - LEVEL 3 (TMR CONTROL 3) 肿瘤标记物质控品 (水平3)						
批号 (Lot. No): 409TU		货号 (Cat. No): TU5003	到期 (Expiry): 2027-04-28					
包装规格 (Size): 3 x 2 ml				Range	范围			
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	ISD	2SD	Method
Alpha-fetoprotein	甲胎蛋白 (AFP)	ng/ml	89.8	71.8	108	9.10	18.2	BioMerieux VIDAS
		KIU/l = IU/ml	74.3	59.4	89.2	7.45	14.9	
		ng/ml	112	89.6	134	11.0	22.0	Roche Cobas e801
		KIU/l = IU/ml	92.7	74.2	111	9.15	18.3	
Beta-2-Microglobulin	β-2-微球蛋白	μg/ml = mg/l	3.02	2.42	3.62	0.300	0.600	BioMerieux VIDAS
CA 125	糖类抗原 125 (CA125)	U/ml	158	126	190	16.0	32.0	BioMerieux VIDAS
		U/ml	131	105	157	13.0	26.0	Roche Cobas e801
CA 15-3	糖类抗原 15-3 (CA15-3)	U/ml	74.9	59.9	89.9	7.50	15.0	BioMerieux VIDAS
		U/ml	87.5	70.0	105	8.75	17.5	Roche Cobas e801
CA 19-9	糖类抗原 19-9 (CA19-9)	U/ml	43.0	34.4	51.6	4.30	8.60	BioMerieux VIDAS
		U/ml	38.2	30.6	45.8	3.80	7.60	Roche Cobas e801
CEA	癌胚抗原	ng/ml = μg/l	40.7	32.6	48.8	4.05	8.10	BioMerieux VIDAS
		ng/ml = μg/l	31.5	25.2	37.8	3.15	6.30	Roche Cobas e801
CYFRA 21-I	细胞角蛋白片段-21	ng/ml	16.0	12.0	20.0	2.00	4.00	Roche Cobas e801
Ferritin	铁蛋白	ng/ml = μg/l	285	228	342	28.5	57.0	BioMerieux VIDAS
		ng/ml = μg/l	268	214	322	27.0	54.0	Randox Immunoturbidimetric assay 朗道免疫比浊分析法
PSA Free	游离前列腺特异抗原 (fPSA)	ng/ml = μg/l	26.8	20.1	33.5	3.35	6.70	Roche Cobas e801
PSA Total	总前列腺特异抗原 (tPSA)	ng/ml = μg/l	37.2	27.9	46.5	4.65	9.30	BioMerieux VIDAS
		ng/ml = μg/l	37.7	28.3	47.1	4.70	9.40	Roche Cobas e801
Thyroglobulin	甲状腺球蛋白	ng/ml	153	115	191	19.0	38.0	Roche Cobas e801
Total Beta hCG	总 β 人绒毛膜促性腺激素	IU/ml	0.213	0.170	0.256	0.022	0.043	BioMerieux VIDAS
		mIU/ml = IU/l	213	170	256	21.5	43.0	
		IU/ml	0.119	0.095	0.143	0.012	0.024	Roche Cobas e801
		mIU/ml = IU/l	119	95.2	143	12.0	24.0	